

Products



RELATED PRODUCTS:

- Diacell® Diamond anvil cells
- 5MP Colour Camera

RELATED ACCESSORIES:

- Eyepieces WF20X/F.N 20
- LED illumination unit
- Polarized light ring

Premium Stereomicroscope

Stereomicroscope ideal for Diamond Anvil Cells sample preparation. The magnification and illumination are compatible with simple work for DACs including diamond anvils alignment .

- ◆ Greenough stereoscopic binocular 450 optical system
- ◆ Trinocular viewing head, inclined 45° with C-mount adapter
- ◆ Working Distance: 115mm
- ◆ Eyepiece: WF10X/F.N.22
- ◆ Magnification Range: 0.8X~5X
- ◆ Dual focusing knobs
- ◆ Dimmable power LED illumination system for iridescent and transmitted light.
- ◆ Optional polarized light available under request
- ◆ An optional camera is available to be mounted on the top of the microscope instead of the observation eyepiece.
- ◆ The camera can be connected to a PC by a USB cable and enable the observation of the gasket hole, sample and Ruby calibrant.

Technical Specifications:

Interpupillary Adjustment	52 to 76 mm
Diopter Adjustment	±5°
Zoom ratio	1:6.3
Magnification range	0.8X-5X
Objective	1X
Working Distance	115 mm
Eyepiece	WF10X/F.N.22
Field of View	27.5mm – 4.4mm

Specifications subject to change without prior notice.
easyLab and Diacell are registered trademarks of Almax easyLab

www.almax-easyLab.com


Almax easyLab bv
Wagenmakerijstraat 5
8600 Diksmuide
Belgium
Ph: +32 51 55 56 37

Almax easyLab Inc (For US and Canada)
Harvard Square -1, Mifflin Place
Cambridge, MA 02138,
United States of America
Ph: + 1 617 701 7245





Stereomicroscope configurations



Stereomicroscope with camera option

